




MICROSCOPY AND
LABORATORY DIAGNOSTICS

Our global presence

“ We manufacture products for use in histopathology, hematology, cytology, and microbiology with a wide national and international distribution network present in over 70 countries.

BioGnost offers only the best quality products at competitive prices and has one of the widest ranges of medical products in South-Eastern Europe. Our experience and flexibility towards the client makes us a sought-after business partner, which has been recognized on the global market.



Our entrepreneurship started in 1990 with the founding of the family company. Accelerated business growth began in 1998 with production of blood typing serums followed by the Croatian Ministry of Health granting BioGnost the License for the manufacture of medical products. We are proud of one of our biggest successes - relocating to our new commercial and manufacturing building in 2012. The site enabled an enhanced work environment, including a carefully devised warehouse structure, segmented production facilities with increased manufacturing capacity, and additional R&D capabilities.

“ We continue to build on the strong momentum from last year to create high-quality products and provide excellent customer service.


Strategic priorities

Our mission is to create high-quality products and provide excellent services, as well as introduce new or expanded product lines and sales ranges in accordance with market dynamics and demands.

Product and service excellence makes BioGnost a desirable regional and international business partner. The key to success is the ever-expanding network of business partnerships, as well as recognizable social responsibility. In the coming years we plan to further our market expansion, additionally invest in research and development, and to become a competitive manufacturer in the fields of microbiological diagnostics and histopathology.

ISO 9001, ISO 13485, ISO 14001

BioGnost's entire business is in accordance with international ISO 9001, ISO 13485 and ISO 14001 standards. The company implemented ISO standards more than 15 years ago, and is the first Croatian company to implement the ISO 13485 quality system for medical devices. Having fulfilled the requirements of these standards, the company is completely in accordance with national legal regulation and the European Union's Medical Device Directive.



CE, IVD

BioGnost's products bear the CE mark, a certificate that confirms that the products are in accordance with the European Union's standards. They are declared as *in vitro* diagnostic medical devices (IVD) classified as "Others" (according to the Directive 98/79/EC). No control of a Notified Body is required for this class; conduct of the entire company business is in accordance with applicable regulations and standards.

OEM manufacture

Besides products from our own sales ranges, we also offer OEM manufacture for your company. It includes a wide range of services, such as label design, product packaging, instructions for use and certificates of analysis, and products manufactured according to your specifications.

Contents



FIXATION AND PROCESSING 6

Marking systems 6

Formalin fixatives 8

Other fixatives 10

Decalcifying solutions 11

Dehydrating and rehydrating agents 12

Clearing and deparaffinization agents 13

Infiltration, embedding and demolding media 14

MOUNTING AND MICROSCOPY 16

Mounting media 16

Immersion oils 20

Cleaning solutions for microscope lens and slides 21

STANDARD HISTOLOGY STAINING REAGENTS 22

Hematoxylin reagents 22

Differentiation and bluing reagents 24

Eosin reagents 25

CYTOLOGY AND HEMATOLOGY STAINING REAGENTS 26

Papanicolaou (Pap) reagents 26

Rapid and standard hematology reagents and kit 28

Cytochemical staining kits and ancillary reagents **NEW** 32

BACTERIOLOGY AND MYCOLOGY STAINING REAGENTS 34

BioGram kits and reagents 34

TB-Stain kits and reagents 36

Solutions for staining and analysis of fungi 39

SPECIAL STAINING KITS FOR HISTOLOGY AND CYTOLOGY 40

ADDITIONAL STAINING REAGENTS AND SOLUTIONS 51

POWDER DYES AND STAINS 58

ACCESSORIES 62



Formalin fixatives

Fixation is one of the most important steps in tissue processing, preventing tissue degradation and decomposition. Preserving tissue properties is important for detection of morphological and pathological changes in tissue samples.

4% NB Formaldehyde

10% neutral buffered and stabilized formalin solution, pH 7.0.

1000 mL	FNB4-1L
5000 mL	FNB4-5L
10000 mL	FNB4-10L
20000 mL	FNB4-20L



4% Formaldehyde

Stabilized 10% formalin solution.

1000 mL	F4-1L
5000 mL	F4-5L
10000 mL	F4-10L
20000 mL	F4-20L



Decalcifying solutions

Decalcifiers most often contain inorganic or organic acids enabling rapid decalcification (i.e. calcium removal). Milder types of decalcifiers are formulated from chelating agents that are used for decalcification of more sensitive types of tissue.

OsteoFast 1

Light blue fixative and decalcifying solution for bone and hard tissue in histology, based on formaldehyde and hydrochloric acid.

1000 mL OF1-OT-1L
2500 mL OF1-OT-2.5L



OsteoFast 2

Light yellow decalcifying solution for bone and hard tissue in histology, based on formic acid.

1000 mL OF2-OT-1L
2500 mL OF2-OT-2.5L



OsteoSens

Decalcifying solution for sensitive calcified tissues in histology, contains EDTA.

1000 mL OS-OT-1L
2500 mL OS-OT-2.5L



Clearing and deparaffinization agents

Clearing agents are completely miscible with alcohol and paraffin, making them an ideal intermedium in histology and cytology processing. They are used multiple times during tissue processing and section staining; before embedding the tissue in paraffin, for deparaffinization before staining and before covering the stained section using medium for permanent preservation.

BioClear

Xylene-based medium for use in histology and cytology.

1000 mL	BC-1L
2500 mL	BC-2.5L
5000 mL	BC-5L
10000 mL	BC-10L

BioClear New

Xylene substitute, aliphatic hydrocarbon-based. Essentially no odor. For use in histology and cytology.

1000 mL	BCN-1L
2500 mL	BCN-2.5L
5000 mL	BCN-5L
10000 mL	BCN-10L



CryoFix Gel

Embedding matrix for cryosectioning of tissue samples known as OCT medium (optimal cutting temperature). Efficient at temperatures of -10°C and below. Available in 5 different colors.

CryoFix Gel	100 mL	CF-100
CryoFix Gel pink	100 mL	CFR-100
CryoFix Gel blue	100 mL	CFP-100
CryoFix Gel green	100 mL	CFZ-100
CryoFix Gel yellow	100 mL	CFY-100



Cryo Spray

Medium in spray bottle for freezing histological samples.

200 mL 08-SPRAY-2



DeMold Spray

Solution in spray bottle for easy removal of paraffin tissue blocks from the mold.

100 mL DS-100



DeMold solution

Solution for refilling DeMold Spray bottle.

500 mL DS-500
1000 mL DS-1L



Tissue Soft

Solution for softening the tissue embedded in paraffin block before microtome sectioning.

500 mL TS-500
1000 mL TS-1L



BioGnost's synthetic BioMount media properties

NAME OF BIOMOUNT MEDIUM	SOLVENT	POLYMER	VISCOSITY @ 20°C (cSt, mPa)	REFRACTIVE INDEX @ 20°C (n _p)	USAGE IN AUTOMATED COVER SLIP MACHINES
BIOMOUNT DPX	xylene (BioClear)	S	850-950	1.5210-1.5230	-
BIOMOUNT DPX HIGH	xylene (BioClear)	S	1900-2100	1.5235-1.5255	-
BIOMOUNT DPX LOW	xylene (BioClear)	S	140-160	1.5115-1.5135	✓
BIOMOUNT DPX NEW*	xylene (BioClear)	S	200-300	1.5150-1.5180	✓
BIOMOUNT NEW	isoparaffin (BioClear New)	M	950-1050	1.4385-1.4395	-
BIOMOUNT NEW LOW	isoparaffin (BioClear New)	M	130-160	1.4345-1.4365	✓
BIOMOUNT	xylene (BioClear)	M	450-650	1.4920-1.4930	✓
BIOMOUNT HIGH	xylene (BioClear)	M	5000-5500	1.4900-1.4920	-
BIOMOUNT M	toluene	M	140-160	1.4900-1.4920	✓
BIOMOUNT C	xylene (BioClear)	M	5000-5500	1.4910-1.4925	-

Legend: M – methacrylate, S – styrene (polystyrene)
 * does not contain harmful dibutyl phthalate plasticizer
 ** contact us to find out how to decrease viscosity of a medium

Xylene and styrene-based media for permanent mounting

BioMount DPX

Moderate viscosity medium for permanent covering of microscope samples and mounting cover slips.

100 mL BM-100
 250 mL BM-250
 500 mL BM-500



BioMount DPX High

High viscosity medium for permanent covering of microscope samples and mounting cover slips.

100 mL BMH-100
 250 mL BMH-250
 500 mL BMH-500



BioMount DPX Low

Low viscosity medium for permanent covering of microscope samples and mounting cover slips. Can be used in automatic cover slip machines.

100 mL BML-100
 250 mL BML-250
 500 mL BML-500



BioMount DPX New

Low viscosity medium for permanent covering of microscope samples and mounting cover slips. Medium contains more environmentally friendly plasticizer instead of harmful dibutyl phthalate. Can be used in automatic cover slip machines.

100 mL BMLE-100
 250 mL BMLE-250
 500 mL BMLE-500





Methacrylate and xylene substitute-based media for permanent mounting

BioMount New

Moderate viscosity medium for permanent covering of microscope samples deparaffinized using xylene substitute-based medium, as well as xylene-based medium.

100 mL	BMN-100
250 mL	BMN-250
500 mL	BMN-500

BioMount New Low

Very low viscosity medium for permanent covering of microscope samples deparaffinized using xylene substitute, as well as xylene.

100 mL	BMNL-100
250 mL	BMNL-250
500 mL	BMNL-500



Methacrylate-based permanent mounting media

BioMount

Low viscosity xylene-based medium for permanent covering of microscope slides and mounting cover slips.

100 mL	BMT-100
250 mL	BMT-250
500 mL	BMT-500

BioMount High

Very high viscosity xylene-based medium for permanent covering of microscope samples and mounting cover slips.

100 mL	BMTH-100
250 mL	BMTH-250
500 mL	BMTH-500

BioMount M

Very low viscosity toluene-based medium for permanent covering of microscope samples and mounting cover slips.

100 mL	BMM-100
250 mL	BMM-250
500 mL	BMM-500



Additional permanent mounting media

BioMount C

Very high viscosity medium for permanent covering of microscope slides and mounting cover slips.

100 mL	BMC-100
250 mL	BMC-250
500 mL	BMC-500

* excellent synthetic substitution for natural covering media, such as Canada Balsam



Water-based mounting media

NAME OF BIOMOUNT MEDIUM	SOLVENT	POLYMER	VISCOSITY @ 20°C (cSt, mPa)	REFRACTIVE INDEX @ 20°C (n _D)	USAGE IN AUTOMATED COVER SLIP MACHINES
BIOMOUNT AQUA	water	PVP	1800-2100	1.3630-1.3670	-



BioMount Aqua

Water-based medium for covering microscope slides and mounting cover slips. With application in immunohistochemistry as a barrier between chromogen stains and the organic solvent-based permanent mounting media.

100 mL	BMA-100
250 mL	BMA-250
500 mL	BMA-500



Immersion oils

Seven different types of immersion oils for microscopical analysis

Immersion oils are essential for routine laboratory work, especially in hematology, cytology and microbiology. They are commonly used for high magnifications in order to secure high resolution and a clear image. There are several types of immersion oils that differ according to viscosity and certain properties required by different microscopy techniques. They all have a common refractive index of 1.5150, equal to the refractive index of glass.

BioGnost's Immersion oils properties

NAME OF IMMERSION OIL	VISCOSITY @ 23°C (cSt, mPa)	COMBINING WITH OTHER TYPES OF OIL	TYPE OF MICROSCOPY
IMMERSION OIL	800-1300	Immersion oil type A Immersion oil type C	light
IMMERSION OIL TYPE A	100-300	Immersion oil Immersion oil type C	light
CEDARWOOD OIL	1500-2200	-	light
IMMERSION OIL TYPE C	1500-2200	Immersion oil Immersion oil type A	light
IMMERSION OIL TYPE FF	170	-	fluorescence
IMMERSION OIL TYPE 37*	1250**	-	light, for use at 37°C
IMMERSION OIL, TROPICAL GRADE***	1000-1500****	-	light, for use at 27-29°C

* refractive index n_D at 37°C = 1.5150

** viscosity at 37°C

*** refractive index n_D at tropical temperatures = 1.515

**** viscosity at 28°C

Immersion oil

Immersion oil for microscopy analysis with optimal viscosity of 600-1100 cSt.

10x10 mL	IU-1010
30 mL	IU-30
100 mL	IU-100
250 mL	IU-250
500 mL	IU-500

Immersion oil type A

Immersion oil for microscopy analysis with low viscosity of 100-300 cSt.

10x10 mL	IUA-1010
30 mL	IUA-30
100 mL	IUA-100
250 mL	IUA-250
500 mL	IUA-500





Immersion oil type C

Immersion oil for microscopy analysis with high viscosity of 1500-2200 cSt.

100 mL	IUC-100
250 mL	IUC-250
500 mL	IUC-500

* excellent synthetic substitution for natural immersion oils, such as Cedarwood oil

Immersion oil type FF

Immersion oil for sample analysis in fluorescence microscopy.

100 mL	IUF-100
250 mL	IUF-250
500 mL	IUF-500

Immersion oil type 37

Immersion oil for microscopy analysis with optimal refractive index and viscosity of 1250 cSt at 37°C.

100 mL	IU37-100
250 mL	IU37-250
500 mL	IU37-500

Immersion oil, tropical grade

Immersion oil for microscopy analysis in tropical climates with optimal refractive index at 27-29°C.

10x10 mL	IUT-1010
30 mL	IUT-30
100 mL	IUT-100
250 mL	IUT-250
500 mL	IUT-500



Cleaning solutions for microscope lens and slides

Clean Lens, solution 1

Solution according to Carl Zeiss formulation for cleaning microscope lens. Based on n-hexane and 2-propanol.

100 mL	CL1-OT-100
--------	------------

Clean Lens, solution 3

Solution for cleaning microscope lens and removing immersion oil from stained slide after using microscope. Based on petroleum ether and 2-propanol.

100 mL	CL3-OT-100
--------	------------

Clean Lens, solution 2

Solution for cleaning microscope lens and removing immersion oil from stained slide after using microscope. Based on diethyl ether and 2-propanol.

100 mL	CL2-OT-100
--------	------------



Standard histology staining reagents


Reagents for hematoxylin and eosin staining, the gold standard for demonstration of cellular and tissue structures, providing detailed staining and contrasting of frozen and paraffin tissue sections, enabling clear visualization of nuclear, cytoplasmic and other cell structures.

Hematoxylin reagents

Nuclear stains for standard HE staining

Properties of BioGnost's Hematoxylin reagents

TYPE OF HEMATOXYLIN	CYTOLOGY	HISTOLOGY	IHC	PROGRESSIVE METHOD	REGRESSIVE METHOD
HEMATOXYLIN H	-	✓	-	✓	✓
HEMATOXYLIN M	-	✓	✓	✓	-
HEMATOXYLIN ML	-	✓	✓	✓	✓
HEMATOXYLIN G1	✓	✓	✓	✓	-
HEMATOXYLIN G2	✓	✓	✓	✓	-
HEMATOXYLIN G3	✓	✓	-	✓	-
HEMATOXYLIN CAR	-	✓	✓	✓	-
HEMATOXYLIN COL	-	✓	✓	✓	✓
HEMATOXYLIN HP, PAP 1A	✓	-	-	✓	✓
HEMATOXYLIN HP, PAP 1B	✓	-	-	✓	✓
HEMATOXYLIN W KIT	-	✓	✓	✓	acid resistant
HEMATOXYLIN P.T.A.	-	✓	-	✓	-



Each of BioGnost's Hematoxylin is specially stabilized in order to maintain the reagent's reactivity for prolonged periods of time, as well as to maintain the quality of the reagent. Every staining provides a crystal clear and highly detailed image of the cell's nucleus, which in turn makes abnormal or pathological changes easier to detect.

Hematoxylin H

Modified hematoxylin according to Harris, most commonly used formulation in HE staining. Strong intensity reagent for progressive and regressive staining in histopathology.

100 mL	HEMH-OT-100
500 mL	HEMH-OT-500
1000 mL	HEMH-OT-1L
2500 mL	HEMH-OT-2.5L

Hematoxylin M

Modified hematoxylin according to Mayer. Weaker intensity reagent for progressive staining in histopathology and counterstaining in immunohistochemistry.

100 mL	HEMM-OT-100
500 mL	HEMM-OT-500
1000 mL	HEMM-OT-1L
2500 mL	HEMM-OT-2.5L

Hematoxylin ML

Modified hematoxylin according to Mayer-Lillie. Strong intensity new generation reagent for progressive and regressive staining in histopathology.

100 mL	HEMML-OT-100
500 mL	HEMML-OT-500
1000 mL	HEMML-OT-1L
2500 mL	HEMML-OT-2.5L

Hematoxylin G1

Hematoxylin according to Gill, ideal for staining goblet cells. New generation reagent for delicate progressive staining in histopathology, cytology and for counterstaining in immunohistochemistry.

100 mL	HEMG1-OT-100
500 mL	HEMG1-OT-500
1000 mL	HEMG1-OT-1L
2500 mL	HEMG1-OT-2.5L

Hematoxylin G2

Hematoxylin according to Gill, contains double concentration of hematoxylin compared to Hematoxylin G1. Moderate intensity new generation reagent for progressive staining in histopathology, cytology and for counterstaining in immunohistochemistry.

100 mL	HEMG2-OT-100
500 mL	HEMG2-OT-500
1000 mL	HEMG2-OT-1L
2500 mL	HEMG2-OT-2.5L

Hematoxylin G3

Modified hematoxylin contains triple concentration of hematoxylin compared to Hematoxylin G1. Strong intensity new generation reagent for progressive staining in histopathology and cytology.

100 mL	HEMG3-OT-100
500 mL	HEMG3-OT-500
1000 mL	HEMG3-OT-1L
2500 mL	HEMG3-OT-2.5L





Hematoxylin reagents

Nuclear stains for standard HE staining

Hematoxylin W kit

Acid-resistant hematoxylin according to Weigert. Two-reagent kit that stains the nuclei blue-black, often a component of special staining kits for connective tissues.

2x100 mL	HEMW-K-100
2x500 mL	HEMW-K-500
2x1000 mL	HEMW-K-1L

Hematoxylin P.T.A.

Synonyms: Hematoxylin Phosphotungstic Acid, Mallory P.T.A.H. reagent. Simultaneously stains different types of tissue reddish-brown and blue. For differentiation between fibrin, smooth and cross-striated muscle fibers.

500 mL	HPTA-OT-500
1000 mL	HPTA-OT-1L

Hematoxylin CAR

Modified hematoxylin according to Carazzi for delicate nuclear staining without staining cytoplasmic components. Reagent for progressive staining of frozen sections in histopathology.

100 mL	HEMCA-OT-100
500 mL	HEMCA-OT-500
1000 mL	HEMCA-OT-1L

Hematoxylin COL

Modified hematoxylin according to Cole, contains iodine. Weaker intensity reagent for progressive and regressive staining in histopathology.

100 mL	HEMCO-OT-100
500 mL	HEMCO-OT-500
1000 mL	HEMCO-OT-1L

Hematoxylin Instant kit

Two-component kit that contains ready-to-use powders for preparation of hematoxylin reagent for progressive and regressive nuclear staining in histopathology and cytology.

for 1000 mL	HEMI-K-1L
for 6000 mL	HEMI-K-6L

Kit for rapid HE staining of frozen and paraffin-embedded sections

Ready-to-use eight-reagent kit (in 20 containers that can be used as staining jars) for rapid HE staining of frozen and paraffin tissue sections in histopathology. Contains xylene substitute as clearing agent and xylene substitute-based medium for permanent coverslipping.

for 100 tests	HE-RTU-100T
---------------	-------------

NEW

ErioGnost reagent

Acid fast reagent for nuclear staining, comparable to hematoxylin staining. Ecologically acceptable synthetic replacement for hematoxylin reagents.

100 mL	ERS-OT-100
500 mL	ERS-OT-500
1000 mL	ERS-OT-1L

Differentiation and bluing reagents for histology and cytology

Differentiation of the nucleus enables clear visualization of nuclear structures in regressive procedure, while the bluing reagent turns the nucleus blue.

Acid alcohol

Reagent for enabling crisp visualization of nuclear chromatin detail and removing of excess background dye in histological samples after regressive staining with hematoxylin.

1000 mL	KA-OT-1L
2500 mL	KA-OT-2.5L

HCL Pap reagent

Buffered reagent for enabling crisp visualization of nuclear chromatin detail and removing of excess background dye in cytological after regressive staining with hematoxylin.

1000 mL	HCLPAP-OT-1L
2500 mL	HCLPAP-OT-2.5L

Scott's solution

Buffered reagent for conversion of red hematoxylin stained nuclei to optimal blue.

1000 mL	SC-OT-1L
2500 mL	SC-OT-2.5L

Bluing reagent

Buffered reagent for conversion of red hematoxylin stained nuclei to optimal blue.

1000 mL	BR-OT-1L
2500 mL	BR-OT-2.5L

Eosin reagents

Cytoplasmic stains for standard HE staining

Eosin Y (yellowish eosin) solutions are commonly used as a counterstain to hematoxylin in the standard histology HE staining method. Eosin Y is an anionic dye that stains basic cell components (such as cytoplasm, collagen, muscle fibers and erythrocytes) bright red.

Eosin Y aqueous solutions

Eosin Y aqueous solutions for cytoplasmic counterstaining.

Eosin Y 0.2% aqueous	1000 mL 2500 mL	EOY-02-OT-1L EOY-02-OT-2.5L
Eosin Y 0.5% aqueous	1000 mL 2500 mL	EOY-05-OT-1L EOY-05-OT-2.5L
Eosin Y 1% aqueous	1000 mL 2500 mL	EOY-10-OT-1L EOY-10-OT-2.5L
Eosin Y 2% aqueous	1000 mL 2500 mL	EOY-20-OT-1L EOY-20-OT-2.5L

Eosin Y alcoholic solutions

Eosin Y alcoholic solutions for cytoplasmic counterstaining.

Eosin Y 0.5% alcoholic	1000 mL 2500 mL	EOYA-05-OT-1L EOYA-05-OT-2.5L
Eosin Y 1% alcoholic	1000 mL 2500 mL	EOYA-10-OT-1L EOYA-10-OT-2.5L
Eosin Y 2% alcoholic	1000 mL 2500 mL	EOYA-20-OT-1L EOYA-20-OT-2.5L

Eosin Contrast

Modified alcoholic solution for cytoplasmic counterstaining. The reagent contains Eosin Y and Phloxine B for enhanced counterstaining effect.

1000 mL 2500 mL	EOYK-OT-1L EOYK-OT-2.5L
--------------------	----------------------------

Eosin Contrast PLUS

Modified alcoholic solution for intensive cytoplasmic counterstaining. Reagent contains Eosin Y, Phloxine B, and Biebrich Scarlet dyes for additional counterstaining effect.

1000 mL 2500 mL	EOYKP-OT-1L EOYKP-OT-2.5L
--------------------	------------------------------

Eosin Y instant

Mixture of powders for preparation of eosin reagent in aqueous or alcoholic medium.

Eosin Y Instant, aqueous	for 1000 mL	EOYI-AQ-1L
Eosin Y Instant, alcoholic	for 1000 mL	EOYI-ALC-1L





Cytology and hematology staining reagents

Cytology and hematology staining kits contain reagents for fixation, staining and buffering cytology and hematology samples, as well as cytology smears of hematopoietic organs, providing precise distinction between cellular structures. Cytochemical staining kits contain reagents for cytochemical diagnosis of leukemia using bone marrow or whole blood smears, and also for staining and determining specific cell components in blood samples.

Papanicolaou (Pap) reagents

Nuclear and cytoplasmic reagents for Pap staining

Because of their simplicity of use and exceptional precision in diagnostics, staining reagents for staining according to Papanicolaou have substantially contributed to decreasing occurrence of cervical cancer. By using a simple staining protocol, BioGnost's Pap reagents enable excellent differentiation between cellular structures.



Hematoxylin HP, Pap 1A

Modified Harris hematoxylin solution for nuclear staining according to Papanicolaou. Strong intensity reagent for progressive and regressive staining in exfoliative cytology.

100 mL	HEMHP-OT-100
500 mL	HEMHP-OT-500
1000 mL	HEMHP-OT-1L
2500 mL	HEMHP-OT-2.5L

Hematoxylin HP, Pap 1B

Hematoxylin modified according to Harris for intensive nuclear staining according to Papanicolaou. Reagent for progressive and regressive staining in exfoliative cytology, with formulation stronger than Hematoxylin HP, Pap 1A.

100 mL	HEMHP1B-OT-100
500 mL	HEMHP1B-OT-500
1000 mL	HEMHP1B-OT-1L
2500 mL	HEMHP1B-OT-2.5L

OG-6 reagent, Pap 2A

Monochromatic reagent for counterstaining mature and keratinized squamous cell. Contains Orange G dye that stains the cytoplasm orange.

100 mL	OG6-OT-100
500 mL	OG6-OT-500
1000 mL	OG6-OT-1L
2500 mL	OG6-OT-2.5L

Orange II reagent, Pap 2B

Monochromatic reagent for counterstaining mature and keratinized squamous cells. Contains Orange II dye that stains the cytoplasm intense red.

100 mL	OR2-OT-100
500 mL	OR2-OT-500
1000 mL	OR2-OT-1L
2500 mL	OR2-OT-2.5L



Rapid and standard hematology reagents and kits

Bio-Diff kits

Rapid staining solutions for hematology, cytology and histology

Rapid diagnostic kits that provide precise staining results in a short period of time, such as results provided by May-Gruenwald Giemsa (MGG) staining method. Each component of Bio-Diff and Bio-Diff RTU kits is stabilized and prepared according to the highest standards.

Bio-Diff kit

Three-reagent kit that contains fixative agent, red and blue components for fast and effective staining. Each kit contains buffer tablets for consistent staining results.

3x100 mL	BD-K-100
3x500 mL	BD-K-500
3x1000 mL	BD-K-1L
3x2500 mL	BD-K-2.5L



Bio-Diff 1 reagent

Fixative agent, component of Bio-Diff kit.

100 mL	BD1-OT-100
500 mL	BD1-OT-500
1000 mL	BD1-OT-1L
2500 mL	BD1-OT-2.5L



Bio-Diff 2 reagent

Buffered red component of Bio-Diff kit.

100 mL	BD2-OT-100
500 mL	BD2-OT-500
1000 mL	BD2-OT-1L
2500 mL	BD2-OT-2.5L

Bio-Diff 3 reagent

Buffered blue component of Bio-Diff kit.

100 mL	BD3-OT-100
500 mL	BD3-OT-500
1000 mL	BD3-OT-1L
2500 mL	BD3-OT-2.5L

STAINING DONE IN 14 SECONDS*



* depending on the type of sample

Bio-Diff RTU kit

Ready-to-Use three-reagent kit with reagents in containers that can be used as staining jars. Kit contains fixative agent, red and blue components for fast and effective staining and buffer tablet for consistent staining results.

3x100 mL BD-RTU-100



Bio-Diff 1 RTU reagent

Fixative component Bio-Diff RTU kit.

100 mL BD1-RTU



Bio-Diff 2 RTU reagent

Buffered red component of Bio-Diff RTU kit.

100 mL BD2-RTU

Bio-Diff 3 RTU reagent

Buffered blue component of Bio-Diff RTU kit.

100 mL BD3-RTU

STAINING DONE IN
14 SECONDS*



* depending on the type of sample

Buffer solution, pH 7.2

Phosphate buffer solution, molarity 6 mM for use in hematology and histology.

1000 mL PUF72-OT-1L
2500 mL PUF72-OT-2.5L
5000 mL PUF72-OT-5L

Buffer tablets pH 7.2

Phosphate buffer tablets for preparing buffer solutions used in hematology and histology. One buffer tablet is sufficient for preparation of 1000 mL of phosphate buffer solution.

50 pcs PT-72-50
100 pcs PT-72-100



Use buffer solution pH 6.8 for hematological and cytological diagnostics. Use buffer solution pH 7.2 for the detection of blood parasites.



May-Gruenwald and Giemsa reagents

Romanowsky staining solutions for hematology, cytology and histology

One of the most commonly used Pappenheim methods is May-Gruenwald Giemsa (MGG) for staining hematology and cytology materials, as well as cytology smears of hematopoietic organs. Giemsa solution is also used for staining histology sections.

May-Gruenwald solution

Polychromatic fixative solution containing Eosin Y, Methylene Blue and azure dyes. For staining in hematology, cytology and staining sections of hematopoietic organs in histopathology.

100 mL	MG-OT-100
10x100 mL	MG-OT-110
500 mL	MG-OT-500
1000 mL	MG-OT-1L
2500 mL	MG-OT-2.5L



Giemsa solution

Polychromatic solution of Eosin Y, Methylene Blue and azure dyes. For use in hematology, cytology, histology and for staining section samples of hematopoietic organs.

100 mL	GM-OT-100
10x100 mL	GM-OT-110
500 mL	GM-OT-500
1000 mL	GM-OT-1L
2500 mL	GM-OT-2.5L



Buffer solution, pH 6.4

Phosphate buffer solution, molarity 6 mM for use in hematology and histology.

1000 mL	PUF64-OT-1L
2500 mL	PUF64-OT-2.5L
5000 mL	PUF64-OT-5L

Buffer solution, pH 6.8

Phosphate buffer solution, molarity 6 mM for use in hematology and histology.

1000 mL	PUF68-OT-1L
2500 mL	PUF68-OT-2.5L
5000 mL	PUF68-OT-5L

Buffer tablets pH 6.8

Phosphate buffer tablets for preparing solutions used in hematology and histology. One buffer tablet is sufficient for preparation of 1000 mL of phosphate buffer solution.

50 pcs	PT-68-50
100 pcs	PT-68-100

Use buffer solution pH 6.8 for hematological and cytological diagnostics. Use buffer solution pH 7.2 for the detection of blood parasites.



Additional staining solutions for hematology and cytology



Wright's solution

Polychromatic solution of Eosin Y, Methylene Blue and azure dyes. For staining in hematology, cytology and cytogenetics.

100 mL	WR-OT-100
500 mL	WR-OT-500
1000 mL	WR-OT-1L

Field kit

Ready to use two-reagent kit for rapid and efficient staining and detection of parasites in hematology samples. Kit for 100 tests contains reagents in containers that can be used as staining jars.

for 100 tests	FD-100T
2x500 mL	FD-K-500

Wright-Giemsa solution

Polychromatic solution of Eosin Y, Methylene Blue and azure dyes. Modified Romanowsky stain for hematology staining.

100 mL	WRGM-OT-100
500 mL	WRGM-OT-500
1000 mL	WRGM-OT-1L
2500 mL	WRGM-OT-2.5L

Hayem's solution

Solution for manual counting of erythrocytes.

100 mL	HY-OT-100
500 mL	HY-OT-500
1000 mL	HY-OT-1L

Tuerk's solution

Solution for manual counting of leukocytes.

100 mL	TU-OT-100
500 mL	TU-OT-500
1000 mL	TU-OT-1L



Brilliant Cresyl Blue, solution

Solution for staining reticulocytes in hematology.

100 mL	BCB-OT-100
500 mL	BCB-OT-500
1000 mL	BCB-OT-1L

Leishman's solution

Polychromatic solution of Eosin Y, Methylene Blue and azure dyes. For staining in hematology and clinical cytology.

100 mL	LE-OT-100
500 mL	LE-OT-500
1000 mL	LE-OT-1L

NEW

CytoGnost kit

Three-reagent kit for manual counting of human blood components

3 x 100 mL	HCG-K-100
3 x 500 mL	HCG-K-500
3 x 1000 mL	HCG-K-1L

ThromboGnost solution

Solution for manual counting of thrombocytes.

100 mL	THR-OT-100
500 mL	THR-OT-500
1000 mL	THR-OT-1L

BioGram ECO kit PHENOL-FREE

Four-reagent phenol-free kit for the identification of bacteria according to Gram. Kit contains Gram Crystal violet, phenol-free reagent, Gram Sodium hydrogencarbonate, solution, stabilized Gram Lugol solution, double amount of Gram Decolorizer solution 2 and Gram Safranin solution as counterstain.

2x50 mL+4x100 mL	BGRE-K-100
2x250 mL+4x500 mL	BGRE-K-500

Gram Crystal Violet, phenol-free reagent

Primary stain reagent for Gram-positive bacteria, doesn't contain phenol. Crystal Violet is also known as Gentian Violet and Methyl Violet 10B.

100 mL	GCF-OT-100
250 mL	GCF-OT-250
500 mL	GCF-OT-500
1000 mL	GCF-OT-1L
2500 mL	GCF-OT-2.5L

Gram Sodium hydrogencarbonate, solution

Solution of sodium hydrogencarbonate for use with BioGram ECO kit.

100 mL	GNHK-OT-100
250 mL	GNHK-OT-250
500 mL	GNHK-OT-500
1000 mL	GNHK-OT-1L
2500 mL	GNHK-OT-2.5L

Additional bacteriology staining solutions according to Gram

Gram Lugol solution

Non-stabilized Lugol solution.

100 mL	GL-OT-100
250 mL	GL-OT-250
500 mL	GL-OT-500
1000 mL	GL-OT-1L

Gram Decolorizer solution 1

Slow acting decolorizer solution used after non-stabilized iodine solution.

100 mL	GD1-OT-100
250 mL	GD1-OT-250
500 mL	GD1-OT-500
1000 mL	GD1-OT-1L

Gram Decolorizer solution 3

Fast acting solution used with automated stainers.

100 mL	GD3-OT-100
250 mL	GD3-OT-250
500 mL	GD3-OT-500
1000 mL	GD3-OT-1L

Gram Safranin 5x concentrated

Concentrated counterstain solution of Safranin O dye.

For preparation of 500 mL of Safranin O counterstain solution	100 mL	SO05X-OT-100
For preparation of 2500 mL of Safranin O counterstain solution	500 mL	SO05X-OT-500
For preparation of 5000 mL of Safranin O counterstain solution	1000 mL	SO05X-OT-1L

Carbol Fuchsin Gram 10x concentrated

Concentrated Carbol Fuchsin solution for preparation of counterstain solution for modified staining. Also known as concentrated Gram Fuchsin solution.

For preparation of 1000 mL of Carbol Fuchsin Gram counterstain solution	100 mL	CFGR10-OT-100
For preparation of 5000 mL of Carbol Fuchsin Gram counterstain solution	500 mL	CFGR10-OT-500
For preparation of 10000 mL of Carbol Fuchsin Gram counterstain solution	1000 mL	CFGR10-OT-1L

Additional iodine solutions

Lugol's solution (Schiller's test)

Classic formulation of aqueous solution of iodine and potassium iodide. Commonly used in gynecology.

100 mL	LUG-OT-100
250 mL	LUG-OT-250
500 mL	LUG-OT-500

Lugol's C solution

Aqueous solution of iodine and potassium iodide, enhanced formulation with stronger staining (compared to classic Lugol solution).

100 mL	LUGC-OT-100
250 mL	LUGC-OT-250
500 mL	LUGC-OT-500

Lugol's solution, acid

Aqueous iodine and potassium iodide solution with added acid.

100 mL	LUGK-OT-100
250 mL	LUGK-OT-250
500 mL	LUGK-OT-500



TB-Stain kits and reagents

Kits for staining acid-fast bacteria according to Ziehl-Neelsen

Most commonly known method for detecting TB bacteria according to Ziehl-Neelsen. This method enables rapid and easy staining of acid-fast bacterial strains, such as *Mycobacteria*, *Nocardia*, *Rhodococcus*, *Dietzia* and *Gordonia*. TB-Stain ECO kit is a phenol-free, environment-friendly alternative to TB-Stain Hot kit.

TB-Stain Hot kit

Three-reagent kit for staining acid-fast bacteria. Contains TB Carbol Fuchsin reagent, double amount of TB Decolorizer and Methylene Blue Loeffler's reagent as counterstain.

for 100 tests	TBH-100T
4x100 mL	TBH-K-100
4x250 mL	TBH-K-250
4x500 mL	TBH-K-500



TB Carbol Fuchsin reagent

Synonym: Ziehl-Neelsen reagent. Primary stain reagent for acid-fast bacteria used in TB-Stain Hot, TB-Stain Cold and TB-Stain Quick kits.

100 mL	TBC-OT-100
250 mL	TBC-OT-250
500 mL	TBC-OT-500
1000 mL	TBC-OT-1L
2500 mL	TBC-OT-2.5L



TB-Stain ECO kit PHENOL-FREE

Three-reagent phenol-free kit for staining acid-fast bacteria. Contains TB-Fuchsin reagent, double amount of TB Decolorizer and Methylene Blue Loeffler's reagent as counterstain.

5x100 mL	TBE-K-100
5x250 mL	TBE-K-250
5x500 mL	TBE-K-500



TB Decolorizer

For use in TB-Stain Hot and TB-Stain Cold kits.

100 mL	TBD-OT-100
250 mL	TBD-OT-250
500 mL	TBD-OT-500
1000 mL	TBD-OT-1L
2500 mL	TBD-OT-2.5L



Methylene Blue Loeffler's reagent

Blue counterstain for TB-Stain Hot kit.

100 mL	MBL-OT-100
250 mL	MBL-OT-250
500 mL	MBL-OT-500
1000 mL	MBL-OT-1L
2500 mL	MBL-OT-2.5L



TB-Fuchsin reagent

Primary stain phenol-free reagent for acid-fast bacteria used in TB-Stain ECO kit.

100 mL	TBFR-OT-100
250 mL	TBFR-OT-250
500 mL	TBFR-OT-500
1000 mL	TBFR-OT-1L
2500 mL	TBFR-OT-2.5L



Additional TB-Stain Hot kit staining solution

Methylene Blue 10x concentrated

Concentrated counterstain solution of Methylene Blue dye for use in microbiology.

For preparation of 1000 mL of Methylene Blue counterstain solution	100 mL	MB10X-OT-100
--	--------	--------------

For preparation of 5000 mL of Methylene Blue counterstain solution	500 mL	MB10X-OT-500
--	--------	--------------

For preparation of 10000 mL of Methylene Blue counterstain solution	1000 mL	MB10X-OT-1L
---	---------	-------------





Kits for detection of acid-fast bacteria using fluorescence microscopy according to modified Truant method

TB-Stain Fluorescent kit

Three-reagent kit for detection of acid-fast bacteria using fluorescence microscopy. Contains TB Auramine-Rhodamine reagent, double amount of TB Decolorizer Fluorescent and TB Permanganate reagent as counterstain.

for 100 tests	TBF-100T
4x100 mL	TBF-K-100
4x250 mL	TBF-K-250
4x500 mL	TBF-K-500



TB-Stain Auramine O kit

Three-reagent kit for staining acid-fast bacteria using fluorescence method. Contains TB Auramine O reagent, double amount of TB Decolorizer Fluorescent and counterstain of TB Permanganate reagent.

4x100 mL	TBAO-K-100
4x250 mL	TBAO-K-250
4x500 mL	TBAO-K-500



TB Auramine O reagent

Primary stain for use in TB Stain Auramine O kit.

100 mL	TBAO-OT-100
250 mL	TBAO-OT-250
500 mL	TBAO-OT-500
1000 mL	TBAO-OT-1L
2500 mL	TBAO-OT-2.5L

TB Auramine-Rhodamine reagent

Primary stain reagent for use in TB-Stain Fluorescent kit.

100 mL	TBA-OT-100
250 mL	TBA-OT-250
500 mL	TBA-OT-500
1000 mL	TBA-OT-1L
2500 mL	TBA-OT-2.5L



TB Decolorizer Fluorescent

For use in TB-Stain Fluorescent and TB-Stain Auramine O kits.

100 mL	TBF-OT-100
250 mL	TBF-OT-250
500 mL	TBF-OT-500
1000 mL	TBF-OT-1L
2500 mL	TBF-OT-2.5L




TB Permanganate reagent

Counterstain for use in TB-Stain Fluorescent kit and TB-Stain Auramine O kits.

100 mL	TBP-OT-100
250 mL	TBP-OT-250
500 mL	TBP-OT-500
1000 mL	TBP-OT-1L
2500 mL	TBP-OT-2.5L






A.F.O.G. kit

Six-reagent Acid Euchsin Orange G kit for selective staining of glomerular protein deposits and collagen in kidney biopsies. Nuclear stain is obtained with Weigert ferric hematoxylin, cytoplasm with Orange G and highly selective collagen stain with Aniline blue.


for 100 tests AFOG-100T
6x100 mL AFOG-K-100



Alcian Blue - P.A.S. kit

Seven-reagent Alcian Blue - Periodic Acid-Schiff kit for staining acid mucopolysaccharides according to Mowry. Enables differentiation between acid mucins (stained light blue) and neutral mucins, glycogens and glycoproteins (stained magenta).


for 100 tests AB-100T
7x100 mL AB-K-100
7x500 mL AB-K-500



Alcian Blue pH 1.0 kit

Three-reagent kit for specific staining of heavily sulfated mucopolysaccharides blue. Slides are counterstained with Nuclear Fast Red (Kernechtrot) reagent to fully detect the presence of Alcian Blue positive staining.


for 100 tests AB10-100T
3x100 mL AB10-K-100



Alcian Blue pH 2.5 kit

Three-reagent kit for staining sulfated and carboxylated acid mucopolysaccharides and sulfated and carboxylated sialomucins blue. Slides are counterstained with Nuclear Fast Red (Kernechtrot) reagent to fully detect the presence of Alcian Blue positive staining.


for 100 tests AB25-100T
3x100 mL AB25-K-100



Alcian Yellow Toluidine Blue kit

Six-reagent kit for staining *Helicobacter pylori* in gastric tissue sections. This method is one of the most popular non-silver method for staining of *H. pylori*, where bacteria are stained blue in contrast to yellow mucins.

for 100 tests AYTB-100T




Azan Trichrome kit

Five-reagent kit for connective tissue staining according to Mallory. Used for visualization of muscle fibers, collagen, glial cells, kidney glomerular stroma, erythrocytes, reticulum and nuclear chromatin in histologic sections.

for 100 tests AZT-100T
5x250 mL AZT-K-250






Elastica-Van Gieson kit

Four-reagent kit for staining elastic fibers and differentiation between elastic tissue, collagen and other types of connective tissue. The rapid method enables a satisfactory result with shorter section staining time.


for 100 tests WVGB-100T
4x100 mL WVGB-K-100



Eosin and nigrosine vital kit

Two-reagent kit containing separate dyes for rapid detection of sperm vitality and simple visualization of dead and living sperm cells. The Nigrosin stain provides dark background for easier recognition of both viable and non-viable spermatozoa.


2x30 mL ENVK-K-30
2x100 mL ENVK-K-100



Feulgen kit

Five-reagent DNA staining kit according to Feulgen. For use in semiquantitative DNA determination in histological and cytological samples. Specimen is first treated with hydrochloric acid creating an aldehyde group of DNA that can be visualized by Schiff (BioSchiff) reagent. This reaction is specific for nuclear DNA.

for 100 tests FE-100T
5x100 mL FE-K-100



Field kit

Ready-to-use two-reagent kit for rapid and efficient staining and detection of parasites in hematology samples. Primarily used for staining thin and thick blood smears (dense drop) for purpose of diagnosing blood parasites. Reagents are stored in containers that can be used as staining jars.


for 100 tests FD-100T



Fouchet-Van Gieson kit

Three-reagent kit for visualization of bilirubin and collagen according to Kutlick. Bilirubin is a yellow-brown pigment, but changes to green due to oxidation induced by Fouchet reaction. Green bilirubin can easily be detected on yellowish and pink colored background.

for 100 tests FVG-100T
3x100 mL FVG-K-100




Giemsa HP kit

Four-reagent kit for staining *Helicobacter pylori* in gastroscopic sections according to Lennart. Advantages of this method for detecting *H. pylori* are sensitive and reproducible results and easy performance.

for 100 tests GMHP-100T






Gomori Trichrome kit

Five-reagent kit for staining muscle, collagen fibre and nuclei, contains blue counterstain. Kit can be used to contrast skeletal, cardiac or smooth muscle.

for 100 tests GT-100T
5x100 mL GT-K-100


* Gomori Trichrome kit with green counterstain also available



Grimelius kit

Five-reagent kit for staining argyrophilic granules. Grimelius kit can be used for the detection of secretory intracytoplasmic granules specific for carcinoid tumors and for identification of neuroendocrine cells.

for 100 tests GRI-100T
5x100 mL GRI-K-100




Grocott kit

Six-reagent kit for visualization of fungi and histological argentaffin structures in general (such as basal membranes). Green counterstain provides clear and visually rich contrast to target structures stained black.

for 100 tests GRC-100T


NEW improved formulation and stability at room temperature!



H.B.F.P. kit

Three-reagent Hematoxylin-Basic Euchsin-Picric acid staining kit for detection of cardiac muscle changes after ischemia or myocardial infarction. H.B.F.P. kit is a non-enzymatic histochemical technique used for identification of early myocardial ischemia occurring on the cardiac muscle during the first six hours of symptoms.


for 100 tests HBFP-100T
3x100 mL HBFP-K-100



Hematoxylin P.T.A. kit

Four-reagent Hematoxylin-Phosphotungstic Acid kit for differentiation of smooth and striated muscle tissues as well as for detection of fibrin, collagen and elements of central nervous system according to Mallory.

for 100 tests HPTA-100T
4x100 mL HPTA-K-100




Hematoxylin W kit

Acid-resistant hematoxylin according to Weigert. Two-reagent kit that stains the nuclei blue-black, often a component of special staining kits for connective tissues.

for 100 tests HEMW-100T
2x100 mL HEMW-K-100
2x500 mL HEMW-K-500
2x1000 mL HEMW-K-1L






HemoGnost Perls kit

Three-reagent HemoGnost Perls (Prussian blue / Berlin blue) kit for the detection of reactive ferric (Fe³⁺) ions in cells. It is often applied on bone marrow and spleen cells.


for 100 tests HP-100T
3x250 mL HP-K-250



Ladewig Trichrome kit NEW

Four-reagent kit for staining muscle (brown to red) and collagen fibers with blue counterstain. Staining procedure is easier and faster compared to Masson and Mallory Trichrome staining procedures.

for 100 tests LDW-100T
4x100 mL LDW-K-100



Luxol Fast Blue kit

Three-reagent kit for staining myelin and myelinated axons, Nissl bodies and phospholipids according to Kluver-Barrera. Kit is used for identifying the basic neuronal structure in brain or spinal cord sections.


for 100 tests LFB-100T
3x100 mL LFB-K-100



M.I.F. kit

Two-reagent Merthiolate-Iodine-Formalin kit for fixing and staining fecal parasites (especially protozoa, cysts, helminth eggs and larvae). Stool sample is fixed with formalin and stained with two coloring agents; iodine and Eosin Y.


for 100 tests MIF-100T



Mallory Trichrome kit

Three-reagent staining kit for connective tissue visualization and detection of collagen, cartilage, muscle, elastic fibers, mucous, pituitary cells, reticulum, bones, amyloid and erythrocytes.

for 100 tests MT-100T
3x100 mL MT-K-100



Martius Scarlet Blue (MSB) kit

Seven-reagent kit used for fibrin visualization, especially of older clusters. This method is a modification of Masson Trichrome and is ideal for studying connective tissue and vascular pathology.

for 100 tests MSB-100T
6x100 mL + 1x250 mL MSB-K-100






Masson Fontana kit

Six-reagent melanine and argentaffin granule staining kit, based on the reduction of silver nitrate to elemental silver. Melanin is a brown-black pigment normally present in hair, skin, retina, iris and certain parts of CNS. Argentaffin granules are found in carcinoid tumors.


for 100 tests MF-100T
9x100 mL MF-K-100



Masson-Goldner Trichrome kit

Seven-reagent kit for staining muscle and collagen fibers with green counterstain. It is also used for visualizing gametes, nuclei, neurofibrils, glial cells, keratins, intercellular fibrils and for differentiation of smooth muscle fibers and collagens.


for 100 tests MGT-100T
7x100 mL MGT-K-100
7x500 mL MGT-K-500



Masson Trichrome kit

Seven-reagent kit for staining muscle and collagen fibers with blue counterstain. It is also used for visualizing gametes, nuclei, neurofibrils, glial cells, keratins and intercellular fibrils. Kit may be useful for differentiation between collagen and smooth muscle fibers or detection of diseases like cirrhosis.


for 100 tests MST-100T
7x100 mL MST-K-100
7x500 mL MST-K-500



Movat kit

Ten-reagent kit for staining collagen, elastic and muscle fibers, mucin and fibrin in tissue section. Movat kit is particularly useful when examining heart and vascular diseases.


for 100 tests MOV-100T
10x100 mL MOV-K-100



Mucicarmine kit

Mucicarmine kit is often used to identify primary tumor sites, distinguishing mucin-negative undifferentiated squamous cell lesions from mucin-positive adenocarcinomas. It can also be used as indicative of diseases such as asthma, bronchitis, chronic obstructive pulmonary disease and cystic fibrosis.

for 100 tests MUC-100T
4x100 mL MUC-K-100




Oil Red O kit

Four-reagent kit for selective staining and detection of fat cells and neutral fats according to Johnson. It can be used with frozen sections and fresh smears to detect obesity-linked pathologies such as dyslipidemia and diabetes.

for 100 tests ORO-100T
4x250 mL ORO-K-250






Orcein kit

Five-reagent kit for visualization of hepatitis B surface antigen (HBsAg) seen as viral inclusion bodies in hepatocytes, for elastic fibers and copper associated protein in tissue sections. It can be used with frozen sections.


for 100 tests ORC-100T
5x250 mL ORC-K-250



Paraldehyde Fuchsin kit

Seven-reagent kit according to Gomori for detecting pathological changes in elastic fibers. It also stains mast cell granules, beta granules in pancreatic islets, neurosecretory material, mast cell granules and beta cells in the pituitary gland.


for 100 tests PAF-100T
9x100 mL PAF-K-100



P.A.S. kit

Five-reagent Periodic Acid-Schiff kit for staining aldehydes, mucopolysaccharides, mucoproteins and lymphocytes according to Hotchkiss-McManus. P.A.S. staining may also be used for the demonstration of fungal organisms in tissue sections.


for 100 tests PAS5-100T
5x100 mL PAS5-K-100
5x500 mL PAS5-K-500



P.A.S. Diastase kit NEW

Six-reagent kit for detecting glycogen and mucopolysaccharide structures mostly in liver. The diastase digestion procedure, followed by P.A.S. staining, is useful in the evaluation of glycogen storage disease.

for 100 tests PDIA-100T




P.A.S.M. / Jones kit

Six-reagent Periodic Acid Silver Methenamine kit for staining kidney glomerular basement membranes. Kit includes red counterstain which provides visually rich contrast to target structures stained black.

for 100 tests PASM-100T
5x100 mL + 3x250 mL PASM-K-100

NEW improved formulation and stability at room temperature!




Perls-Van Gieson kit NEW

Three-reagent kit for the detection of reactive ferric (Fe^{3+}) ions in cells (Prussian blue/Berlin blue reaction) with Fuchsin Acid Van Gieson counterstain to differentiate collagen from connective tissue

for 100 tests PEG-100T






Picro-Sirius Red kit NEW

Three-reagent kit for collagen staining and visualization in histology. Contains Sirius Red dye which enhances collagen birefringence so it can easily be identified by observing the section under polarized light.


for 100 tests PSR-100T



Reticulin kit

Seven-reagent kit for detecting argyrophilic reticulin fibers. It clearly differentiates between collagen and reticulin and nerve fibers and connective tissue. The main function of reticular fibers is to provide support and are normally found in liver, lymph nodes, spleen and kidneys.


for 100 tests RE-100T
7x50 mL RE-K-50
7x100 mL RE-K-100



Reticulin Contrast kit

Nine-reagent kit for detecting argyrophilic reticulin fibers according to Gordon and Sweets. The kit contains gold chloride solution that enhances visualization of reticulin fibers and it also contains Nuclear Fast Red (Kernechtrot) reagent that enables fine contrasting background.


for 100 tests RET-100T
9x50 mL RET-K-50
9x100 mL RET-K-100



Rhodanine kit

Four-reagent kit for detecting copper and copper-associated protein (CAP) in patients suffering from Wilson's disease. Abnormal copper accumulations are predominately found in liver tissue, but can also be found in the brain and cornea.


for 100 tests ROD-100T
5x100 mL ROD-K-100



Safranin O kit NEW

Four-reagent kit for identification of cartilage, mucins and mastocyte granules used in histology and cytology.

for 100 tests SAFO-100T
4x100 mL SAFO-K-100




Sperm-Diff RTU kit

Ready-to-Use three-reagent kit with reagents stored in containers that can be used as staining jars. Kit contains fixative reagent, red and blue components for fast and effective staining.

3x100 mL SP-RTU-100






Sudan Black B Lipid kit

Four-reagent kit for specific lipid staining in cytochemistry. Sudan Black dye stains a few types of lipids, including neutral fats, phospholipids and sterols. Contains double amount of Sudan Black B, solution.


for 100 tests SBB-100T
6x100 mL SBB-K-100



TB-Stain Histo kit

Three-reagent kit for staining acid-fast bacteria (pathogenic mycobacteria) in histology sections, sputum, smears and culture smears according to Ziehl-Neelsen. Heating of the carbol-fuchsin solution is avoided in this protocol hence the release of hazardous phenolic vapors is omitted.

for 100 tests TBHK-100T



UriGnost S kit

Kit for sampling, staining, and microscopic analysis of urine sediment. Contains modified reagent according to Sternheimer (reagent manufactured according to the European Confederation of Laboratory Medicine (ECLM) and all the necessary equipment for sampling, concentrating, staining, and counting cells and kidney cylinders, as well as urine sediment analysis (test tubes with retentive bottom, Eppendorf 200 µL pipette tips and urine plates).


for 500 tests USK-500
* for more information, see page 57



UriGnost SM kit

Kit for sampling, staining and microscopic analysis of urine sediments. Contains modified reagent according to Sternheimer-Malbin and equipment required for sampling, concentration, staining, counting of kidney cells and kidney cylinders and urine sediment analysis (test tubes with retentive bottom, Eppendorf 200 µL pipette tips and urine plates).


for 500 tests USMK-500
* for more information, see page 57



Van Gieson Trichrome kit

Three-reagent kit for staining collagen fibers, muscle tissue, keratinized epithelium, cytoplasm, glial fibers and erythrocytes. Used for differentiation between collagen and smooth fibers in tumors and various other diseases.

for 100 tests VG-100T
3x100 mL VG-K-100




Verhoeff kit

Six-reagent kit for detecting atrophy of elastic tissue in cases of emphysema, thinning and loss of elastic fibers in arteriosclerosis and other vascular diseases, or whether blood vessels have been invaded by a tumor.

for 100 tests VER-100T
6x100 mL VER-K-100






Von Kossa kit

Five-reagent kit for simple and reproducible detection of calcium deposits and calcium salt in tissue samples according to Von Kossa. Tissue calcification is associated with metabolic problems in various tissues (like bone marrow and mamma) and in tumors.


for 100 tests VK-100T
 5x100 mL VK-K-100



Warthin Starry kit

Five-reagent kit for staining *Spirochaeta*, *Helicobacter pylori*, *Microsporadia* and *Legionella pneumophila*. The kit contains 12 jars with gelatin that enables both incubation and staining of sections, as well as other reagents that enable precipitation of silver on the bacterial surface. The bacteria are found in the mucus of the surface epithelium, in the apical gastric glands and in the gastric mucosa.

up to 48 tests WS-12T



Weigert-Van Gieson kit

Six-reagent kit for staining elastic fibers with prolonged incubation period. Used for differentiation between elastic tissue, collagen and other types of connective tissue.

for 100 tests WVG-100T
 6x100 mL WVG-K-100



Additional staining reagents and solutions

Solutions and reagents (buffers, special stains, developer solutions and reagents for differentiation or counterstaining) that are either standalone reagents or components of special staining kits.

Staining reagents and solutions

4% Formaldehyde, alcoholic solution

For use with Reticulin and Reticulin Contrast staining kits.

100 mL F4A-OT-100



Acetic buffer, Grimelius

For use with Grimelius kit

100 mL APG-OT-100

Acid alcohols

For use with Azan Trichrome and Orcein staining kits.

Acid alcohol, Azan ● 100 mL KAA-OT-100
250 mL KAA-OT-250

Acid alcohol, Orcein ● 250 mL KAO-OT-250

Acid solution

For preparation of impregnation solution with Warthin Starry kit.

500 mL KOT-OT-500

Acridine Orange solution

Acridine Orange aqueous solution is used as fluorescent dye for differential staining of DNA and RNA, acid mucins and detecting apoptosis; in cytogenetics for displaying structures containing abundant amounts of DNA. Also enables detection of microorganisms in blood smears.

100 mL AOA-OT-100
500 mL AOA-OT-500
1000 mL AOA-OT-1L

A.F.O.G. reagent

For use with A.F.O.G. (Acid Fuchsin Orange G) kit.

100 mL AFOG-OT-100

Alcian Blue solution 1%

Alcian Blue 8GX dye 1% aqueous solution.

100 mL AB1-OT-100
500 mL AB1-OT-500

Alcian Blue solution pH 1.0

Reagent for staining heavily sulfated mucopolysaccharides. A component of Alcian Blue pH 1.0 kit.

100 mL AB10-OT-100
500 mL AB10-OT-500
1000 mL AB10-OT-1L

Alcian Blue solution pH 2.5

Reagent for staining acid mucopolysaccharides of intestinal goblet cells. A component of Alcian Blue-P.A.S. and Alcian Blue pH 2.5 staining kits.

100 mL AB2-OT-100
500 mL AB2-OT-500
1000 mL AB2-OT-1L

Alkaline alcohol, solution

For use with MSB kit.

100 mL ALA-OT-100



Ammonium water NEW

For use with Bielschowsky kit.

100 mL AV-OT-100

Ammonium iron sulfate, solution

For use with Reticulin and Reticulin Contrast staining kits.

100 mL ASF-OT-100



Aniline, alcoholic solution

For use with Azan Trichrome kit.

100 mL	ANA-OT-100
250 mL	ANA-OT-250

Aniline Blue reagent

Counterstain used with Masson Trichrome kit.

100 mL	ABR-OT-100
500 mL	ABR-OT-500

Aniline Blue, solution

For use with MSB kit.

100 mL	ABO-OT-100
--------	------------

Azan reagent

For use with Azan Trichrome kit.

100 mL	AZR-OT-100
250 mL	AZR-OT-250

Azocarmine, solution

For use with Azan Trichrome kit.

100 mL	AZ-OT-100
250 mL	AZ-OT-250

Basic Activation buffer

For use with Oil Red O kit.

250 mL	BAP-OT-250
--------	------------

Benedict's reagent

For detecting glucose in urine.

100 mL	BNR-OT-100
250 mL	BNR-OT-250
500 mL	BNR-OT-500
1000 mL	BNR-OT-1L

Biebrich Scarlet-Acid Fuchsin reagent

For use with special kits.

100 mL	BSAF-OT-100
500 mL	BSAF-OT-500

BioSchiff reagent

Schiff's reagent for detecting aldehydes and mucous substances. A component of P.A.S. and Alcian Blue-P.A.S. staining kits.

100 mL	BS-OT-100
500 mL	BS-OT-500
1000 mL	BS-OT-1L

Borax, solution

For use with Grocott and P.A.S.M./Jones kits.

100 mL	BO-OT-100
250 mL	BO-OT-250

Chromotrope 2R / Aniline Blue reagent

Blue counterstain used with Gomori Trichrome kit.

100 mL	CABR-OT-100
--------	-------------

Chromotrope 2R / Fast Green F.C.F. reagent

Green counterstain used with Gomori Trichrome kit.

100 mL	CFG-OT-100
--------	------------

Colloidal Iron, stock solution

Reagent for binding ferric ions to acid mucopolysaccharides and proteoglycans.

100 mL	KOL-OT-100
--------	------------

Congo Red Puchtler reagent

For use with Congo Red Puchtler kit.

100 mL	CRPR-OT-100
--------	-------------

SCHIFF'S REAGENT
STABLE AT 22°C



Congo Red reagent

For use with Congo Red Highman kit.

100 mL CR-OT-100



Cresyl Violet, solution

For use with Luxol Fast Blue kit.

100 mL CV-OT-100

Developer stock solution NEW

For use with Bielschowsky kit.

50 mL STR-OT-50

Differentiation reagents based on acetic acid

For use with special staining kits in histopathology.

Reagent for differentiation in Gomori Trichrome kit	100 mL	RDG-OT-100
Acetic acid, 0.5% solution	100 mL	OKO5-OT-100
Acetic acid, 1% solution	100 mL 500 mL	OK1-OT-100 OK1-OT-500
Acetic acid, Colloidal Iron (A) solution	250 mL	OKA-OT-250
Acetic acid, Colloidal Iron (B) solution	100 mL	OKB-OT-100

Differentiation reagents based on hydrochloric acid

For use with special staining kits in histopathology.

HCL reagent, A.F.O.G. ●	100 mL	HCLAF-OT-100
HCL reagent, Alcian Blue	100 mL	HCLAB-OT-100
HCL reagent, Colloidal Iron ●	100 mL	HCLK-OT-100
HCL reagent, Feulgen ●	100 mL	HCLF-OT-100
HCL reagent, HemoGnost Perls	100 mL	HCLH-OT-100
HCL reagent, HemoGnost Perls	250 mL	HCLH-OT-250
HCL reagent, P.A.S.	100 mL	HCLP-OT-100
HCL reagent, P.A.S.	500 mL	HCLP-OT-500
HCL reagent, Weigert-Van Gieson ●	100 mL	HCLW-OT-100

Ehrlich aldehyde reagent

Reagent for detecting urobilinogen in urine.

100 mL EAR-OT-100

Eosin-nigrosin vital

Fast detection (one-step detection) of sperm vitality and visualization of dead and living sperm cells with one reagent. Simple, easy and fast method for semen analysis.

30 mL ENV-30
100 mL ENV-100

Eosin Y, isotonic solution

Aqueous solution of Eosin Y used with Eosin and nigrosin vital kit.

100 mL EIO-OT-100

ErioGnost reagent NEW

Acid fast reagent for nuclear staining, comparable to hematoxylin staining. Ecologically acceptable synthetic replacement for hematoxylin reagents.

100 mL ERS-OT-100
500 mL ERS-OT-500
1000 mL ERS-OT-1L

Fast Green F.C.F. reagent

Counterstain used with Masson-Goldner Trichrome and Paraldehyde Fuchsin staining kits.

100 mL FGR-OT-100
500 mL FGR-OT-500

Fast Green F.C.F. contrast reagent

Counterstain used with Grocott staining kit.

100 mL FGKR-OT-100

Field, solution A

Blue component of Field kit.

100 mL FDA-RTU

Field, solution B

Red component of Field kit.

100 mL FDB-RTU

Ferri reagent, Verhoeff B

For use with Verhoeff kit.

100 mL FRV-OT-100

Ferri reagent, Weigert B

Iron(III) nitrate solution in diluted hydrochloric acid, a component of Hematoxylin W kit.

100 mL FR-OT-100
500 mL FR-OT-500
1000 mL FR-OT-1L

Fuchsin Acid reagent

For use with Mallory Trichrome kit.

100 mL FA-OT-100

Fuchsin Acid Van Gieson reagent

For staining with Van Gieson Trichrome, Weigert-Van Gieson and Colloidal Iron kits. Synonym: Picrofuchsin reagent.

100 mL FAG-OT-100

Fuchsin Basic reagent

For use in bacteriology.

100 mL FBR-OT-100
500 mL FBR-OT-500
1000 mL FBR-OT-1L



Fuchsin Basic, solution

For use with H.B.F.P. kit.

100 mL FBO-OT-100

Gelatin, gel in transport jar

For preparing developing solution with Warthin Starry staining kit.

12x8 mL GEL-128



Glutaraldehyde, 4% buffered solution

Solution used with Sudan Black B kit.

100 mL G4-OT-100





Gold chloride solutions

Aqueous solutions of gold chloride (aurichloric acid, HAuCl_4), component of many special staining kits.

Gold chloride, 0.2% solution	100 mL	ZK02-OT-100
Gold chloride, 0.6% solution	100 mL	ZK06-OT-100
Gold chloride, 1% solution	50 mL	ZK1-OT-50
Gold chloride, 2% solution	50 mL	ZK2-OT-50

Hematoxylin, Verhoeff A

For use with Verhoeff kit.

100 mL	HEMV-OT-100
--------	-------------



Hematoxylin, Weigert A

Hematoxylin solution in ethanol, component of Hematoxylin W kit.

100 mL	HEMA-OT-100
500 mL	HEMA-OT-500
1000 mL	HEMA-OT-1L



Hydroquinone, solution

For preparation of developing solution in Warthin Starry kit.

50 mL	HQ-OT-50
-------	----------



Hydroquinone-sodium sulphite NEW

For use with Grimelius kit

15 x 3,5 g	HNS-15-35
------------	-----------

Iodine solution, Verhoeff C

For use with Verhoeff kit.

100 mL	JODV-OT-100
--------	-------------

Iron chloride, solution

For use with Fouchet-Van Gieson kit.

100 mL	FK-OT-100
--------	-----------

Iron nitrate, solution

Iron nitrate solution for use with Movat kit.

100 mL	FN-OT-100
--------	-----------

Kovac indole reagent

Reagent for identification of indole-positive and indole-negative microorganisms.

100 mL	KOV-OT-100
--------	------------



Ladewig's reagent

For use with Ladewig Trichrome kit, contains Aniline Blue, Acid Fuchsin and Orange G dyes.

100 mL	LDWR-OT-100
--------	-------------

Light Green S.F. reagent

Counterstain used with trichrome special staining kits.

100 mL	LGS-OT-100
--------	------------

Lithium carbonate, Luxol

Lithium carbonate solution used with Luxol Fast Blue kit.

100 mL	LKL-OT-100
--------	------------

Lithium carbonate, solution

Lithium carbonate solution used with Von-Kossa staining kit.

100 mL	LK-OT-100
--------	-----------

Lugol's solution, M.I.F.

For use with M.I.F. kit.

100 mL	LUGM-OT-100
--------	-------------

Luxol Fast Blue, solution

For use with Luxol Fast Blue kit.

100 mL	LFB-OT-100
--------	------------



M.F. reagent

For use with M.I.F. kit.

100 mL	MF-OT-100
--------	-----------

Martius Yellow, solution

For use with MSB kit.

100 mL	MAY-OT-100
--------	------------

Mayer's mucicarmine

For use with Mucicarmine kit.

100 mL	MMU-OT-100
--------	------------



Metanil Yellow, solution

For use with Mucicarmine kit.

100 mL	MY-OT-100
--------	-----------

Methenamine, solution

For use with Grocott and P.A.S.M./Jones kits.

100 mL	MET-OT-100
250 mL	MET-OT-250

Methyl Blue Eosin, solution

Modified Mann's solution for use in histology.

100 mL	MBE-OT-100
500 mL	MBE-OT-500
1000 mL	MBE-OT-1L

Neutral Red reagent

Neutral Red aqueous counterstain solution for use with special staining kits.

100 mL	NR-OT-100
--------	-----------

Nigrosine, isotonic solution

Aqueous solution of nigrosine used with Eosin and nigrosine vital kit.

100 mL	NIO-OT-100
--------	------------

Nuclear Fast Red (Kernechtrot) reagent

Nuclear Fast Red dye solution for staining nuclei. Counterstain in many special staining kits.

100 mL	KR-OT-100
250 mL	KR-OT-250
500 mL	KR-OT-500

Oil Red O reagent

Lipid substances staining solution.

250 mL	ORO-OT-250
--------	------------



Orange G, 1% solution

For use with Movat kit.

100 mL ORG1-OT-100

Orange G / Aniline Blue reagent

For use with Mallory Trichrome kit.

100 mL OGA-OT-100

Orcein reagent

For use with Orcein kit.

250 mL ORCR-OT-250

Oxalic acid solutions

For use with Hematoxylin P.T.A., Masson Fontana, Orcein, Paraldehyde Fuchsin and Reticulin staining kits.

Oxalic acid, 1% solution 100 mL OKS1-OT-100

Oxalic acid, 2% solution 250 mL OKS2-OT-250

Paraldehyde Fuchsin reagent

Reagent for staining elastins and mucins.

100 mL PAF-OT-100

Periodic acid, solution

Synonym: Orthoperiodic acid. Oxidizing medium widely used in histology. Key component of many special staining kits, such as P.A.S. and Alcian Blue-P.A.S., P.A.S.M. / Jones, Grocott and Weigert-Van Gieson kit.

Periodic acid, 0.5% solution 100 mL PK05-OT-100

Periodic acid, 0.5% solution 250 mL PK05-OT-250

Periodic acid, 0.5% solution 500 mL PK05-OT-500

Periodic acid, 0.8% solution 100 mL PK08-OT-100

Periodic acid, 0.8% solution 250 mL PK08-OT-250

Periodic acid, 0.8% solution 500 mL PK08-OT-500

Periodic acid, 1% solution 100 mL PK1-OT-100

Periodic acid, 1% solution 250 mL PK1-OT-250

Periodic acid, 1% solution 500 mL PK1-OT-500

Phosphate buffered saline (PBS)

Sterilized isotonic phosphate buffer solutions.

Buffered saline (Dulbecco), without added calcium and magnesium, sterile filtered 500 mL D-OT-500

Buffered saline (Dulbecco), without added calcium and magnesium, concentrated 10x, sterile filtered 500 mL DX-OT-500

Buffered saline (Dulbecco), with added calcium and magnesium, sterile filtered 500 mL DCM-OT-500

Buffered saline (Dulbecco), with added calcium and magnesium, concentrated 10x, sterile filtered 500 mL DCMX-OT-500

Buffered saline (Dulbecco) (0.01 mol/L, Tween 20), sterile filtered 500 mL PBS-OT-500

Picro-Sirius Red reagent NEW

For use with Picro-Sirius Red kit for collagen staining and differentiation.

100 mL PSR-OT-100

Phosphomolybdic acid, 1% solution

For use with A.F.O.G. (Acid Fuchsin Orange G) and Mallory Trichrome staining kits.

100 mL FMK1-OT-100

Phosphomolybdic acid, 5% solution

For use with Azan Trichrome kit.

100 mL FMK5-OT-100

250 mL FMK5-OT-250

Phosphotungstic acid, 1% solution

For use with Mallory Trichrome kit.

100 mL FVK1-OT-100

500 mL FVK1-OT-500

Phosphotungstic acid, 5% solution

For use with Movat kit.

100 mL FVK5-OT-100

Picric acid in acetone, solution

For use with H.B.F.P. kit.

100 mL PKA-OT-100

Ponceau S, solution

For use with MSB kit.

100 mL PONS-OT-100

Potassium hexacyanoferrate, Colloidal Iron solution

Potassium hexacyanoferrate (II) solution for use with Colloidal Iron kit.

100 mL KHCK-OT-100

Potassium hexacyanoferrate, solution

Potassium hexacyanoferrate (II) solution, optimally concentrated for use with HemoGnost Perls kit.

100 mL KHC-OT-100

250 mL KHC-OT-250

Potassium hydroxide, solution

Optimally concentrated aqueous solution of potassium hydroxide used with Congo Red Highman kit.

100 mL KHO-OT-100

Potassium permanganate, solutions

For use with Hematoxylin P.T.A., Masson Fontana, Orcein, Paraldehyde Fuchsin and Reticulin staining kits.

Potassium permanganate, 100 mL KP05-OT-100
0.5% solution ●

Potassium permanganate, 250 mL KP1-OT-250
1% solution ●



Propylene glycol

For use with Sudan Black B Lipid kit.

100 mL PG-100

Propylene glycol, 85%

For use with Sudan Black B Lipid kit.

100 mL PG85-100

P.T.A. - P.M.A. reagent

Reagent containing optimal concentration of phosphotungstic and phosphomolybdic acid. Used with Masson Trichrome and Masson-Goldner Trichrome staining kits.

100 mL PPR-OT-100
500 mL PPR-OT-500

Reagent for differentiation in Verhoeff kit

For use with Verhoeff kit.

100 mL RDV-OT-100

Reducing buffer, solution

For use with Von Kossa kit.

100 mL RP-OT-100

Resorcine Fuchsin reagent

For use with Weigert-Van Gieson kit.

100 mL RFR-OT-100

Rhodanine reagent

For use with Rhodanine kit.

100 mL RR-OT-100

Safranin O reagent

Counterstaining reagent for use in histology and cytology. Used for visualizing and staining cartilage, mucins and mastocytes in formalin-fixed and paraffin-embedded tissue samples.

100 mL SAFO-OT-100

Silver nitrate solutions

For use in special staining kits.

Silver ammonia solution (component of Reticulin and Reticulin Contrast staining kits) ●● SA-OT-100

Silver ammonia reagent (component of Von Kossa and Masson Fontana staining kits) ●● SAR-OT-100

Silver ammonia reagent stock solution (component of Bielschowsky staining kit) SSAR-OT-100

Silver nitrate, 20% solution (component of Bielschowsky staining kit) SN20-OT-100

Silver nitrate solution, stabilized (component of Grocott and P.A.S.M. / Jones staining kits) ●● SNS-OT-100

Silver nitrate, WS solution (component of Warthin Starry staining kit) ● SNWS-OT-100

Sodium acetate, solution

For use with Rhodanine kit.

100 mL NA-OT-100

Sodium hydroxide, solution

For use with Congo Red Puchtler kit.

30 mL NH-OT-30

Sodium metabisulfite solutions

Component of Alcian Blue-P.A.S., P.A.S. and Feulgen staining kits.

Sodium metabisulfite 100 mL NM-OT-100
500 mL NM-OT-500

Sodium metabisulfite, solution for Feulgen 100 mL NMF-OT-100

Sodium tetraborate, solution

A component of Alcian Blue pH 2.5 and Alcian Blue-P.A.S. staining kits.

100 mL NTB-OT-100
500 mL NTB-OT-500

Sodium thiosulfate solution

Solutions are component of Grocott, Masson Fontana, P.A.S.M. / Jones, Reticulin, Reticulin Contrast and Von Kossa staining kits.

Sodium thiosulfate, 2% solution 100 mL NT2-OT-100

Sodium thiosulfate, 5% solution 100 mL NT5-OT-100

Sudan Black B, solution

Solution for staining lipids, used with Sudan Black B lipid kit.

100 mL SBO-OT-100

Sudan Black B reagent

Staining reagent for use in hematology, used with Sudan Black B kit.

100 mL SBB-OT-100

Sulfosalicylic acid, 20% solution

Reagent for detecting proteins in urine.

100 mL SSK-OT-100

Sulfuric acid solutions

Used with Orcein, Paraldehyde Fuchsin and Reticulin staining kits.

Sulfuric acid, 0.3% solution 250 mL SK03-OT-250

Sulfuric acid, 0.5% solution 100 mL SK05-OT-100

Sulfuric acid, 3% solution 100 mL SK3-OT-100

Tincture of iodine, 1%

Alcoholic iodine solution.

100 mL TJ-OT-100
1000 mL TJ-OT-1L



Hematoxylin, C.I. 75290

Natural Black 1, BSC certified stain. For preparing nuclear staining reagents.

25 g H-P-25
100 g H-P-100

Janus Green B, C.I. 11050

Diazine Green, BSC certified stain. For displaying mitochondria in living cells.

25 g JG-P-25

Jenner Stain

Jenner, BSC certified stain. Romanowsky polychromatic stain routinely used in hematology for staining blood smears.

25 g JS-P-25

Leishman Stain

Eosin-polychrome Methylene Blue. Romanowsky polychromatic dye. Used to distinguish and identify leukocytes and malaria parasites.

25 g LS-P-25

Light Green S.F., C.I. 42095

Acid Green 5, Pencil Green SF, BSC certified stain. For counterstaining according to Papanicolaou and Masson trichrome methods.

25 g LG-P-25

Luxol Fast Blue MBSN, C.I. 74180

Luxol Fast Blue, Solvent Blue 38. For use in histology.

25 g LFB-P-25

Malachite Green, C.I. 42000

Malachite Green Oxalate, Basic Green 4, Victoria Green, BSC certified stain. Used as counterstain in TB-Stain Cold kit.

10 g MAG-P-10
25 g MAG-P-25

Martius Yellow C.I. 10315

Golden Yellow, Martinsgelb, Acid yellow 24. For staining erythrocytes yellow in trichrome staining methods.

25 g MAY-P-25

May-Gruenwald

Polychromatic Romanowsky dye. For blood smears and bone marrow staining.

25 g MG-P-25
50 g MG-P-50

Metanil Yellow, C.I. 13065

Acid yellow 36. For staining cytoplasm, muscle connectors and identifying parietal cells of gastric glands.

25 g MEY-P-25
50 g MEY-P-50
100 g MEY-P-100

Methyl Green, C.I. 42590

For staining nuclei green.

25 g MGR-P-25

Methyl Violet 2B, C.I. 42535

Basic Violet 1, BSC certified stain. For bacterial staining according to Gram.

25 g MV2B-P-25

Methyl Violet 10B, C.I. 42555

Gentian Violet, Crystal Violet, BSC certified stain. For bacterial staining according to Gram.

25 g MV10B-P-25
50 g MV10B-P-50

Methylene Blue, C.I. 52015

Basic Blue 9, BSC certified stain. Used for preparing Romanowsky polychromatic stain for staining bone marrow and as a counterstain in staining of acid-fast bacteria.

25 g MB-P-25
100 g MB-P-100

Methylene Green, C.I. 52020

Methylene Green, BSC certified stain. For staining nuclei green.

25 g MEG-P-25

Neutral Red, C.I. 50040

Basic Red 5, BSC certified stain. For supravital staining of blood cells using Neutral Red-Janus Green B method.

25 g NR-P-25

Nigrosine WS, C.I. 50420

Nigrosine water soluble, Nigrosine B. For use in microbiology and cytology, as negative staining.

25 g NGR-P-25

Nuclear Fast Red, C.I. 60760

Kernechtrot. Used for staining nuclei red by creating aluminum complex.

5 g NF-P-5

Oil Red O, C.I. 26125

Sudan Red 5B, BSC certified stain. For staining lipids.

25 g ORO-P-25

Orange II, C.I. 15510

Acid Orange 7, BSC certified stain. Used as a counterstain, as a substitute for Orange G.

25 g O2-P-25

Orange G, C.I. 16230

Acid Orange 10, Acid Fast Orange G, BSC certified stain. For staining keratin according to Papanicolaou, counterstain for hematoxylin.

25 g OG-P-25
100 g OG-P-100

Orcein, C.I. (ed.1) 1242

Natural Red 28, Orcein Synthetic. For staining elastic tissues and viral inclusion bodies of the hepatitis B surface antigen.

5 g ORC-P-5
10 g ORC-P-10

Phloxine B, C.I. 45410

Acid Red 92, Eosin 10B, Eosin Blue, BSC certified stain. For Methylene Blue-Phloxine staining and staining keratin and mucus according to Kreyberg.

25 g PB-P-25

Ponceau S, C.I. 27195

Acid Red 112, Fast Ponceau 2B. For use in electrophoresis.

25 g PS-P-25

Potassium permanganate

Synonym: Potassium (VII) manganate, Condy's crystals. Potassium permanganate for use in microscopy.

50 g KPH-P-50
100 g KPH-P-100

Pyronine B, C.I. 45010

Pyronine B, BSC certified stain. For bacterial, mold and RNA staining.

25 g PYB-P-25

Resorcinol, C.I. 76505

Resorcin, 1,3-benzendiol. Used for differentiation between aldoses and ketoses in Seliwanoff's test by staining ketose sugars dark red.

25 g RES-P-25

Rhodamine B, C.I. 45170

Basic Violet, Acid Red 52. For staining acid-fast bacteria.

25 g RHB-P-25

Rhodanine

DMAB-Rhodanine. For staining deposits of copper bound to proteins.

5 g RHO-P-5

Rose Bengal, C.I. 45440

Acid Red 94, BSC certified stain. For staining in microscopy and histology.

25 g RB-P-25

Safranin O, C.I. 50240

Basic Red 2, Safranin T, Phenazinium, BSC certified stain. Counterstain in the Gram method for detecting Gram-negative bacteria.

25 g SO-P-25

Sirius Red, C.I. 35780

Direct Red 80. Used for staining collagen and amyloid in histology.

25 g SR-P-25

Sudan III, C.I. 26100

Solvent Red 23, Sudan G, Fat Ponceau G, BSC certified stain. For lipid staining according to Lillie-Ashburn.

25 g S3-P-25

Sudan IV, C.I. 26105

Solvent Red 24, Scarlet Red, Biebrich Scarlet R, Fat Ponceau. For staining lipids according to Churukian's modification of Lillie-Ashburn's method.

25 g S4-P-25

Sudan Black B, C.I. 26150

Sudan Black, Solvent Black 3, Fat Black HB. For staining lipids and other hydrophobic materials.

5 g SB-P-5
10 g SB-P-10

Toluidine Blue O, C.I. 52040

Basic Blue 17, Tolonium Chloride, BSC certified stain. For staining mast cells, cartilage and Nissl bodies in the spinal cord.

25 g TB-P-25
100 g TB-P-100

Victoria Blue B, C.I. 44045

Basic Blue 26, Fat Blue B. For use in cytology and histology.

25 g VBB-P-25

Victoria Blue R, C.I. 42563

Basic Blue 8, Fat Blue 4R. For use in microbiology and histology.

25 g VBR-P-25

Water Blue, C.I. 42755

Acid Blue 22, BSC certified stain. For staining connective tissue according to Mallory and Gomori trichrome methods.

25 g WB-P-25

Wright's Stain

Colorants de Wright, BSC certified stain. Part of polychromatic Romanowsky stains commonly used in hematology.

25 g WS-P-25

Xylidine Ponceau, C.I. 16150

Ponceau 2R, Ponceau de Xylidine, BSC certified stain. For Masson trichrome staining.

25 g XP-P-25





Accessories


High quality disposable laboratory equipment including items made of glass, plastic and metal. Product range includes VitroGnost slides for microscopy, sample transport containers and various other accessories.

VitroGnost® Plus Ultra adhesive microscope slides

These slides are ideal for pathology, cytology, and special immunohistochemistry (IHC). VitroGnost Plus Ultra adhesive slide is a new generation slide with a denser and stronger positive charge on the surface to achieve even better adhesion of the sample. VitroGnost Plus Ultra adhesive slide can be used for many scientific purposes, especially with tissues difficult to fix, like breast tissue and samples of brain and skin. Dimensions: 26x76x1 mm.

72 pcs

103001A72



Lid for cassettes

Cassette embedding stainless steel lid, standard.

material	packaging	product code
INOX	1 pc	51050103





Bio/Cut microtome blades - high profile

Bio/Cut H35C microtome blades for cryotome, high profile (wide). Ideal for frozen sections. Suitable for both soft and hard tissues.

50 pcs Bio/Cut H35C

Histoblades Duraedge histology blades - high profile

HistoBlades Duraedge blades with Teflon coating - high profile.

50 pcs 7310

Bio/Cut microtome blades - low profile

Bio/Cut microtome blades manufactured in order to provide most precise sectioning.

Bio/Cut 22A microtome blades for ultra-thin sections, kidneys and lymph nodes, narrow blade 50 pcs Bio/Cut 22A

Bio/Cut 35A microtome blades for soft and hard tissues, low profile 50 pcs Bio/Cut 35A

Bio/Cut 35B microtome blades for soft tissues, low profile 50 pcs Bio/Cut 35B

Bio/Cut 35N microtome blades for thin sections, low profile 50 pcs Bio/Cut 35N

Bio/Cut 35R microtome blades for hard tissues, low profile 50 pcs Bio/Cut 35R

Bio/Cut 35S microtome blades for hard tissues, low profile 50 pcs Bio/Cut 35S

Histoblades Duraedge histology blades - low profile

HistoBlades Duraedge blades with Teflon coating - low profile.

50 pcs 7223

Histoblades Encore histology blades - low profile

Extremely sharp new generation Teflon coated blades - low profile.

50 pcs 7380



Contact and ordering information

BioGnost Ltd.,
Medjugorska 59, Zagreb
10040 Croatia (Europe)

WORKING HOURS
Monday through Friday
8:00 to 16:00 (CET)

e-mail: sales@biognost.hr
website: www.biognost.com
phone: +385 1 2409 997
fax: + 385 1 2404 039



Technical support

Technical questions can be answered by emailing our sales representatives at sandra.ban@biognost.hr, pavle.loncar@biognost.hr, viktoria.liber@biognost.hr and valentina.perkovic@biognost.hr.

They will provide you with any MSDS or product information required, also available on our website at www.biognost.com.

Frequently Asked Questions (FAQ) section is available on our website.



Shipping information

Shipping and handling charges are prepaid and added to the invoice. If requested, we can use your account number or carrier of your choice, and also ship on your third-party freight accounts (Incoterms 2020). Charges vary according to destination, weight and content, and are available upon request.

Our employees are trained to classify, package, mark and label products identified as dangerous goods according to ADR, IMDG and IATA to ensure safe transport.



Index

4% Formaldehyde	8
4% Formaldehyde, alcoholic solution	51
4% NB Formaldehyde	8
10% Formaldehyde	9
10% NB Formaldehyde	9

A

Acetic acid, 0.5% solution	53
Acetic acid, 1% solution	53
Acetic acid, Colloidal Iron (A) solution	53
Acetic acid, Colloidal Iron (B) solution	53
Acetic buffer, Grimelius	51
Acid alcohol	24
Acid alcohols	51
Acid solution	51
Acridine Orange, C.I. 46005	58
Acridine Orange solution	51
A.F.O.G. kit	41
A.F.O.G. reagent	51
Alcian Blue 8 GX, C.I. 74240	58
Alcian Blue - P.A.S. kit	41
Alcian Blue pH 1.0 kit	41
Alcian Blue pH 2.5 kit	41
Alcian Blue solution 1%	51
Alcian Blue solution pH 1.0	51
Alcian Blue solution pH 2.5	51
Alcian Green	58
Alcian Yellow, C.I. 12840	58
Alcian Yellow Toluidine Blue kit	41
Alcoholic Eosin Y solutions	25
Alkaline alcohol, solution	51
Amido Black 10B, C.I. 20470	58
Ammonium iron sulfate, solution	51
Ammonium water	51
Aniline, alcoholic solution	52
Aniline Blue, C.I. 42780	58
Aniline Blue reagent	52
Aniline Blue, solution	52
Auramine O, C.I. 41000	58
Azan reagent	52
Azan Trichrome kit	41
Azocarmine G, C.I. 50085	58
Azocarmine, solution	52
Azophloxine, C.I. 18050	58
Azure A, C.I. 52005	58
Azure B, C.I. 52010	58
Azure C, C.I. 52002	58
Azure II	58
Azure II Eosinate	59

B

Basic Activation buffer	52
Benedict's reagent	52
Biebrich Scarlet-Acid Fuchsin reagent	52
Biebrich Scarlet, C.I. 26905	59
Bielschowsky kit	42
BioClear	13
BioClear New	13
Bio/Cut microtome blades - high profile	64
Bio/Cut microtome blades - low profile	64
Bio-Diff 1 reagent	28
Bio-Diff 1 RTU reagent	29
Bio-Diff 2 reagent	28
Bio-Diff 2 RTU reagent	29
Bio-Diff 3 reagent	28
Bio-Diff 3 RTU reagent	29
Bio-Diff kit	28
Bio-Diff RTU kit	29, 42
BioFix	9

BioFix C	9
BioFix GL	9
BioFix HD	9
BioForm S	9
BioGram 4 kit	34
BioGram ECO kit	35
BioGram Histo kit	42
BioMount	18
BioMount Aqua	19, 33
BioMount C	19
BioMount DPX	17
BioMount DPX High	17
BioMount DPX Low	17
BioMount DPX New	17
BioMount High	18
BioMount M	18
BioMount New	18
BioMount New Low	18
Biopsy pads for cassettes, blue	63
BioSchiff reagent	52
BioWax 52/54	14
BioWax 56/58	14
BioWax Blue	14
BioWax Plus 56/58	14
Bismarck Brown Y, C.I. 21000	59
Bluing reagent	24
Borax, solution	52
Bouin's solution	10
Brilliant Cresyl Blue, C.I. 51010	59
Brilliant Cresyl Blue, solution	31
Brilliant Green Certified, C.I. 42040	59
Buffer solution, pH 6.4	30
Buffer solution, pH 6.8	30
Buffer solution, pH 7.2	29
Buffer tablets pH 6.8	30
Buffer tablets pH 7.2	29

C

Carbol Fuchsin Gram 10x concentrated	35
Carmine, C.I. 75470	59
Carminic acid, C.I. 75470	59
Celestine Blue, C.I. 51050	59
Chromotrope 2R / Aniline Blue reagent	52
Chromotrope 2R, C.I. 16570	59
Chromotrope 2R / Fast Green F.C.F. reagent	52
CitoFix	10
CitoSpray	10
CitoSpray solution	10
Clean Lens, solution 1	21
Clean Lens, solution 2	21
Clean Lens, solution 3	21
Colloidal Iron kit	42
Colloidal Iron, stock solution	52
Congo Red, C.I. 22120	59
Congo Red Highman kit	42
Congo Red Puchtler kit	42
Congo Red Puchtler reagent	52
Congo Red reagent	53
Cresyl Violet Acetate Certified	59
Cresyl Violet, solution	53
CryoFix Gel	15
Cryo Spray	15
Crystal Violet, C.I. 42555	59
CytoGnost kit	31

D

DeHematin	9
DeMold solution	15
DeMold Spray	15
Developer stock solution	53
Differentiation reagents based on acetic acid	53
Differentiation reagents based on hydrochloric acid	53
Duoline marker pen	7

E	
EA 31 reagent, Pap 3A	27
EA 50 reagent, Pap 3B	27
EA 65 reagent, Pap 3C	27
EA 65 reagent, Pap 3D	27
Ehrlich aldehyde reagent	53
Elastica-Van Gieson kit	43
Eosin and nigrosine vital kit	43
Eosin B, C.I. 45400	59
Eosin Contrast	25
Eosin Contrast PLUS	25
Eosin-nigrosin vital	53
Eosin Y aqueous solutions	25
Eosin Y, C.I. 45380	59
Eosin Y instant	25
Eosin Y, isotonic solution	53
Eriochrome Cyanine R, C.I. 43820	59
ErioGnost reagent	24, 53
F	
Fast Green F.C.F., C.I. 42053	59
Fast Green F.C.F. contrast reagent	53
Fast Green F.C.F. reagent	53
Ferri reagents, Verhoeff B	53
Ferri reagent, Weigert B	53
Feulgen kit	43
Field kit	31, 43
Field, solution A	53
Field, solution B	53
Formaldehyde 37% for use in histology	9
Fouchet-Van Gieson kit	43
Fuchsin Acid, C.I. 42685	59
Fuchsin Acid reagent	53
Fuchsin Acid Van Gieson reagent	53
Fuchsin Basic, C.I. 42510	59
Fuchsin Basic reagent	53
Fuchsin Basic - Schiff, C.I. 42500	59
Fuchsin Basic, solution	53
G	
Gabbet reagent	37
Gelatin, gel in transport jar	53
Gentian Violet, C.I. 42555	59
Giemsa	59
Giemsa HP kit	43
Giemsa solution	30
Glutaraldehyde, 4% buffered solution	53
Gold chloride solutions	54
Gomori Trichrome kit	44
Gram Crystal Violet 1% solution	34
Gram Crystal Violet, phenol-free reagent	35
Gram Decolorizer solution 1	35
Gram Decolorizer solution 2	34
Gram Decolorizer solution 3	35
Gram Lugol solution	35
Gram Lugol solution, stabilized	34
Gram Safranin 5x concentrated	35
Gram Safranin solution	34
Gram sodium hydrogencarbonate, solution	35
Grimelius kit	44
Grocott kit	44
H	
Hayem's solution	31
H.B.F.P. kit	44
HCL Pap reagent	24
HCL reagent, A.F.O.G.	53
HCL reagent, Alcian Blue	53
HCL reagent, Colloidal Iron	53
HCL reagent, Feulgen	53
HCL reagent, HemoGnost Perls	53
HCL reagent, P.A.S.	53
HCL reagent, Weigert-Van Gieson	53
Hematoxylin CAR	24
Hematoxylin, C.I. 75290	60
Hematoxylin COL	24
Hematoxylin G1	23, 27
Hematoxylin G2	23, 27
Hematoxylin G3	23, 27
Hematoxylin H	23
Hematoxylin HP, Pap 1A	26
Hematoxylin HP, Pap 1B	26
Hematoxylin Instant kit	24
Hematoxylin M	23
Hematoxylin ML	23
Hematoxylin P.T.A.	24
Hematoxylin P.T.A. kit	44
Hematoxylin, Verhoeff A	54
Hematoxylin, Weigert A	54
Hematoxylin W kit	24, 44
HemoGnost Perls kit	45
Histanol 50	12
Histanol 70	12
Histanol 80	12
Histanol 95	12
Histanol 96	12
Histanol 100	12
Histanol EP 95	12
Histanol EP 100	12
Histanol G	12
Histanol IP	12
Histanol M	12
Histoblades Duraedge histology blades - high profile	64
Histoblades Duraedge histology blades - low profile	64
Histoblades Encore histology blades - low profile	64
Hydroquinone-sodium sulphite	54
Hydroquinone, solution	54
I	
Immersion oil	20
Immersion oil, tropical grade	21
Immersion oil type 37	21
Immersion oil type A	20
Immersion oil type C	21
Immersion oil type FF	21
ImmunoForm	9
Iodine solution, Verhoeff C	54
Iron chloride, solution	54
Iron nitrate, solution	54
J	
Janus Green B, C.I. 11050	60
Jenner Stain	60
K	
Kit for rapid HE staining of frozen and paraffin-embedded sections	24
Kovac indole reagent	54
L	
Lactophenol Blue solution	39
Lactophenol solution	39
Ladewig's reagent	54
Ladewig Trichrome kit	45
Leishman's solution	31
Leishman Stain	60
LeukoGnost ACP	32
LeukoGnost ALP	32
LeukoGnost Fixative	33
LeukoGnost Fixative ECO	33
LeukoGnost Fixative GL	33
LeukoGnost HEM	33

LeukoGnost MPO	32	Orange G / Aniline Blue reagent	55
LeukoGnost NSE	32	Orange G, C.I. 16230	60
LeukoGnost PAS	32	Orange II, C.I. 15510	60
LeukoGnost PERLS	32	Orange II reagent, Pap 2B	26
LeukoGnost PLUS	33	Orcein, C.I. (ed.1) 1242	60
LeukoGnost SPE	32	Orcein kit	47
LeukoGnost SPENSE	32	Orcein reagent	55
Lid for cassettes	62	OsteoFast 1	11
Light Green S.F., C.I. 42095	60	OsteoFast 2	11
Light Green S.F. reagent	54	OsteoSens	11
Liquid Blocker marker pen	7	Oxalic acid solutions	55
Lithium carbonate, Luxol	54		
Lithium carbonate, solution	54	P	
Lugol's C solution	35	Paraldehyde Fuchsin kit	47
Lugol's solution, acid	35	Paraldehyde Fuchsin reagent	55
Lugol's solution, M.I.F.	54	P.A.S. Diastase kit	47
Lugol's solution (Schiller's test)	35	P.A.S. kit	47
Luxol Fast Blue kit	45	P.A.S.M. / Jones kit	47
Luxol Fast Blue MBSN, C.I. 74180	60	Periodic acid, solution	55
Luxol Fast Blue, solution	54	Perls-Van Gieson kit	47
		Phloxine B, C.I. 45410	61
M		Phosphate buffered saline (PBS)	55
Malachite Green, C.I. 42000	60	Phosphomolybdic acid, 1% solution	55
Mallory Trichrome kit	45	Phosphomolybdic acid, 5% solution	55
Martius Scarlet Blue (MSB) kit	45	Phosphotungstic acid, 1% solution	55
Martius Yellow C.I. 10315	60	Phosphotungstic acid, 5% solution	55
Martius Yellow, solution	54	Picric acid in acetone, solution	55
Masson Fontana kit	46	Picro-Sirius Red kit	48
Masson-Goldner Trichrome kit	46	Picro-Sirius Red reagent	55
Masson Trichrome kit	46	Polypropylene containers for transporting samples	65
Mayer's mucicarmine	54	Polypropylene dissecting boards	63
May-Gruenwald	60	Ponceau S, C.I. 27195	61
May-Gruenwald solution	30	Ponceau S, solution	55
Metal histology mold	63	Potassium hexacyanoferrate, Colloidal Iron solution	55
Metanil Yellow, C.I. 13065	60	Potassium hexacyanoferrate, solution	55
Metanil Yellow, solution	54	Potassium hydroxide, solution	55
Methenamine, solution	54	Potassium hydroxide, solution for fungi	39
Methyl Blue Eosin, solution	54	Potassium permanganate	61
Methylene Blue 10X concentrated	36	Potassium permanganate, solutions	55
Methylene Blue, C.I. 52015	60	Propylene glycol	56
Methylene Blue Loeffler's reagent	36	Propylene glycol, 85%	56
Methylene Green, C.I. 52020	60	P.T.A. - P.M.A. reagent	56
Methyl Green, C.I. 42590	60	Pyronine B, C.I. 45010	61
Methyl Violet 2B, C.I. 42535	60		
Methyl Violet 10B, C.I. 42555	60	R	
M.F. reagent	54	Rapid Pap staining kit	27
Micro Clean	14	Reagent for differentiation in Gomori Trichrome kit	53
Microtome oil	65	Reagent for differentiation in Verhoeff kit	56
M.I.F. kit	45	Reducing buffer, solution	56
Movat kit	46	Resorcine Fuchsin reagent	56
Mucicarmine kit	46	Resorcinol, C.I. 76505	61
Myco-Mount Blue	39	Reticulin Contrast kit	48
Myco-Mount Red	39	Reticulin kit	48
		Rhodamine B, C.I. 45170	61
N		Rhodanine	61
NB Formaldehyde concentrated 5x	9	Rhodanine kit	48
Neutral Red, C.I. 50040	60	Rhodanine reagent	56
Neutral Red reagent	54	Rose Bengal, C.I. 45440	61
Nigrosine, isotonic solution	54		
Nigrosine WS, C.I. 50420	60	S	
Nuclear Fast Red, C.I. 60760	60	Sacconanno's fixative	10
Nuclear Fast Red (Kernechtrot) reagent	54	Safranin O, C.I. 50240	61
		Safranin O kit	48
O		Safranin O reagent	56
OG-6 reagent, Pap 2A	26	Scott's solution	24
OG-EA31 Pap reagent	27	Shorr's solution	27
OG-EA50 Pap reagent	27	Silver nitrate solutions	56
Oil Red O, C.I. 26125	60	Sirius Red, C.I. 35780	61
Oil Red O kit	46	Sodium acetate, solution	56
Oil Red O reagent	54	Sodium hydroxide, solution	56
Orange G, 1% solution	55	Sodium metabisulfite solutions	56

Sodium tetraborate, solution	56
Sodium thiosulfate solution	56
Sperm-Diff RTU kit	48
Sudan Black B, C.I. 26150	61
Sudan Black B ECO	32
Sudan Black B kit	32
Sudan Black B lipid kit	49
Sudan Black B reagent	56
Sudan Black B, solution	56
Sudan III, C.I. 26100	61
Sudan IV, C.I. 26105	61
Sulfuric acid solutions	56
Sulphosalicylic acid, 20% solution	56

T

TB Armand reagent	37
TB Auramine O reagent	38
TB Auramine-Rhodamine reagent	38
TB Carbol Fuchsin reagent	36, 37
TB Decolorizer	36, 37
TB Decolorizer Fluorescent	38
TB-Fuchsin reagent	36
TB Malachite Green reagent	37
TB Permanganate reagent	38
TB-Stain Auramine O kit	38
TB-Stain Cold kit	37
TB-Stain ECO kit	36
TB-Stain Fluorescent kit	38
TB-Stain Histo kit	49
TB-Stain Hot kit	36
TB-Stain Quick kit	37
ThromboGnost solution	31
Tincture of iodine, 1%	56
Tissue marking dyes	7
Tissue marking dyes, 5 colour set	7
Tissue marking dyes, 7 colour set	7
Tissue Soft	15
Toluidine Blue O, C.I. 52040	61
Toluidine Blue, solution	57
Trichloroacetic acid, solution	57
Tuerk's solution	31

U

UriGnost S kit	49
UriGnost SM kit	49
UriGnost SM reagent	57
UriGnost S reagent	57

V

Van Gieson Trichrome kit	49
Verhoeff kit	49
Victoria Blue B, C.I. 44045	61
Victoria Blue R, C.I. 42563	61
VitroGnost® Plus Ultra adhesive microscope slides	62
Von Kossa kit	50

W

Warthin Starry kit	50
Water Blue, C.I. 42755	61
Weigert-Van Gieson kit	50
Wright-Giemsa solution	31
Wright-Giemsa, solution for reptiles	31
Wright's solution	31
Wright's Stain	61

X

Xylidine Ponceau, C.I. 16150	61
------------------------------	----

Be a part of
our global network
of satisfied customers

BioGnost Ltd.

For development, manufacture and
marketing of pharmaceuticals

Medjugorska 59, 10040 Zagreb, Croatia, EU
phone: +385 1 2409 997, fax: +385 1 2404 039
sales@biognost.hr • www.biognost.com

